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Name ^	Version	Repository	Company
✓ Diligent LINX (Control Arduino, R	3.0.1.192	NI LabVIEW Tools Network	Diligent
Diligent WaveForms VIs	1.0.3.26	NI LabVIEW Tools Network	Diligent inc
Distributed State Machines (DSM)	0.0.0.28	NI LabVIEW Tools Network	Ovak Technologies
DMC GUI Suite	7.0.1.1	NI LabVIEW Tools Network	DMC
DQMH	6.0.0.83	NI LabVIEW Tools Network	DQMH Consortium
DQMH Event Scripter	6.0.0.122	NI LabVIEW Tools Network	DQMH Consortium
DQMH Palette	6.0.0.14	NI LabVIEW Tools Network	DQMH Consortium
DQMH Project Template	6.0.0.91	NI LabVIEW Tools Network	DQMH Consortium
DQMH Thermal Chamber Examp	6.0.0.49	NI LabVIEW Tools Network	DQMH Consortium
DSM Framework Example	1.0.0.5	NI LabVIEW Tools Network	Ovak Technologies
Dual Port SPI Example for LabVIEW	1.0.0.1	NI LabVIEW Tools Network	National Instruments
DVR Helper	1.0.0.14	VIPM Community	EUV Tech
Dynamic Signal Analyzer	1.0.0.35	NI LabVIEW Tools Network	MIT
Dynamic VI Handler	1.5.0.23	NI LabVIEW Tools Network	PhysicsConsult
EAL for LabVIEW	2.0.0.10	NI LabVIEW Tools Network	Bosch Rexroth AG
EasyXML Toolkit for LabVIEW	3.0.0.170	VIPM Community	JKI
Enclose With Error Case Structure	1.0.0.1	VIPM Community	Bienkowski Michal
Encryption Compendium for Lab	2.2.0.13	NI LabVIEW Tools Network	Lvs-Tools.Co.Uk
Enhanced State Diagram Toolkit	1.0.0.7	NI LabVIEW Tools Network	National Instruments
Epoch Date & Time	1.0.0.2	VIPM Community	LabVIEW Open Source Project
Epson Robotics Library	1.1.0.47	NI LabVIEW Tools Network	ImagingLab
Epson Robotics Library	3.2.0.84	NI LabVIEW Tools Network	DigiMetrix
Error Logger	1.0.0.2	NI LabVIEW Tools Network	CPE Systems NZ Ltd
Error Structure Drawer	2.2.0.2	VIPM Community	Stuart J. Box
Eulogy	1.0.1.52	NI LabVIEW Tools Network	Field R&D Services
Event Source Actor	1.1.0.5	NI LabVIEW Tools Network	NI
Expression Scanner	1.0.0.10	NI LabVIEW Tools Network	GROVF
Extended Tiny Encryption Algorit	1.0.0.2	NI LabVIEW Tools Network	National Instruments
Fanuc Robotics Library	1.9.3.2	NI LabVIEW Tools Network	DigiMetrix
Flat UI Controls Suite	1.0.0.1	NI LabVIEW Tools Network	Ovak Technologies
Flat User Interface Controls Suite	1.0.0.1	NI LabVIEW Tools Network	Ovak Technologies
Flatline Controls	1.9.1.33	NI LabVIEW Tools Network	JDP Science
Flexible eHS Real Time Power Elec	2016.0.0.1	NI LabVIEW Tools Network	OPAL-RT Technologies
FPGA IP (IPNet): Digital Buses anc	1.0.0.4	NI LabVIEW Tools Network	National Instruments
FPGA IP (IPNet): LabVIEW FPGA E	1.0.0.4	NI LabVIEW Tools Network	National Instruments
Frame APIs for PDC	1.0.0.2	NI LabVIEW Tools Network	National Instruments
FTDI Driver API for Serial Commu	1.0.3.20	NI LabVIEW Tools Network	Ovak Technologies
FTDI Driver for I2C SPI GPIO	1.0.0.31	NI LabVIEW Tools Network	AI

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# Utilizando a biblioteca LINX

The screenshot shows the LabVIEW Package Manager interface. The search bar contains 'linx'. The results table is as follows:

Name ^	Version	Repository	Company
Digilent LINX (Control Arduino, R	3.0.1.192	NI LabVIEW Tools Network	Digilent
LINX Raspberry Pi Addons	1.0.0.1	VIPM Community	MediaMongrels
NI LabVIEW LINX Toolkit	1.0.0.9	NI LabVIEW Tools Network	NI
Xilinx ChipScopePro Debugging E	2.1.0.11	NI LabVIEW Tools Network	B&A Engineering Inc. (BAI)


The 'Install' button for the selected package is highlighted with a red box.

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



# Utilizando a biblioteca LINX

LINX - Blink (Simple).vi Front Panel

File Edit View Project Operate Tools Window Help

18pt Application Font

### Manual Blink Example

This example demonstrates how to blink an LED on a LINX device by clicking on an LED control in LabVIEW.

**Instructions**

1. Select the **Serial Port** associated with the LINX Device.
2. Select the **Digital Output Channel** connected to the LED.
3. Click the **Run Arrow**.

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